

Cleaning and Decontamination Log, Pre-amp room

At a minimum cleaning will entail the following:

- All counter tops cleaned with a solution of bleach or approved substitute (i.e. Dispatch)
- All extraction hoods cleaned with a solution of approx. 10% bleach or approved substitute (i.e. Dispatch) followed by an ethanol wipe down
- All centrifuges cleaned with a solution of approx. 10% bleach or approved substitute (i.e. Dispatch) followed by an ethanol wipe down
- Replace sticky floor mats at both entry doors to the lab.
- Decontaminate the epMotion deck with Decon Quat, followed by 70% ethanol
- All Qiagubes cleaned with Decon Quat, followed by a water wipe down
- All pipettes cleaned with Decon Quat, followed by an ethanol wipe down (Note: Do not spray pipettes directly, spray Decon Quat or ethanol onto a kimwipe and use the kimwipe to clean the pipettes)
- All door knobs, lab chairs, carts, pens, and telephones, freezer and refrigerator doors/handles, mouse, keypad, cabinet and drawer handles cleaned, with a solution of approx. 10% bleach or approved substitute (i.e. Dispatch)
- Used lab coats to laundry
- Run maintenance procedure on Hamilton STARlet (called weekly maintenance on desktop)

Date	Initials	Date	Initials

Document ID	Revision	Approval	Date Published
1800	7	Supervising Criminalist - Biology	2/18/2019 3:23:16 PM

Cleaning and Decontamination Log, Post-amp rooms

At a minimum cleaning will entail the following:

- All counter tops cleaned with a solution of approx. 10% bleach or approved substitute (i.e. Dispatch)
- All centrifuges cleaned with a solution of approx. 10% bleach or approved substitute (i.e. Dispatch) followed by an ethanol wipe down
- Replace sticky floor mat at entry door to the lab.
- Decontaminate the epMotion deck with Decon Quat, followed by 70% ethanol
- All pipettes cleaned with Decon Quat, followed by an ethanol wipe down (Note: Do not spray pipettes directly, spray Decon Quat or ethanol onto a kimwipe and use the kimwipe to clean the pipettes)
- All door knobs, lab chairs, carts, pens, and telephones, freezer and refrigerator doors/handles, mouse, keypad, cabinet and drawer handles cleaned, with a solution of approx. 10% bleach or approved substitute (i.e. Dispatch)
- Used lab coats to laundry

NOTE: DO NOT USE BLEACH OR ETHANOL ON ANY 3130 PARTS

Date	Initials	Date	Initials

Document ID	Revision	Approval	Date Published
1800	7	Supervising Criminalist - Biology	2/18/2019 3:23:16 PM